

EMDR EQUIPMENT EUROPE

Part of: Albero di Olivo Ltd

PROS AND CONS

IN THE USE OF EMDR EQUIPMENT



Hand/Arm moving or Tapping
to activate stimulate the
Bilateral System



Audio Scan, Tactile Unit or Eye
Scan used to stimulate the
Bilateral System

EMDR EQUIPMENT EUROPE

COMPANY HISTORY

Albero di Olivo Ltd



EMDR Equipment Europe
My Bookstation
Blue Stallion Publication
Plan Bee – Bee Kind

- ▶ Since 2002 we work together with KarlWare in Germany to produce, promote and sell EMDR Equipment
- ▶ 2010 Albero di Olivo Ltd has been found as a New Umbrella
- ▶ EMDR Equipment Europe is one part of the Limited Company
- ▶ In 2014 we started to produce in the UK

EMDR EQUIPMENT EUROPE

HAND/ARM MOVEMENT, TAPPING

- ▶ Instant connection with the Client
- ▶ instant adjustment of speed and intensity
- ▶ Very individual tailored
- ▶ Close contact with the client



- ▶ With the Arm movement Therapist suffer from repetitive strain injury
- ▶ Therapist focus on own action or pain coping rather than on clients reaction
- ▶ Clients can feel intimidated by the hand movement
- ▶ Tapping can be inappropriate to use with clients who have been physical or sexually abused

DOWN SIDE

BILATERAL STIMULATION USING EMDR EQUIPMENT



▶ Audio Scan/Tactile Unit

- ▶ Audio Stimulation and Tactile stimulation are synchronic
- ▶ Therapist can change speed, volume and intensity
- ▶ Easy to control
- ▶ After your client is familiar with the treatment, you can hand over the Control Unit to your Client. This helps on the recovery and is empowering.

BILATERAL STIMULATION USING EMDR EQUIPMENT

▶ Eye Scan Unit

- ▶ Frees your Arm, you can focus more on the processing and how your clients is responding, take notes
- ▶ **You as the Therapist can choose from:**
 - ▶ - different colours red/green/blue/white
 - ▶ - three different settings
 - ▶ - adjust the speed
 - ▶ - control the brightness/dimmable



EMDR EQUIPMENT EUROPE

- ▶ Some Clients find it very technical
- ▶ Some clients maybe will find that there are treated by a machine
- ▶ You have to carry the equipment with you

DOWN SIDE

► **Comprises of:**

The Basic Unit and the Headphones



Basic Unit is the Control Unit

- **Speed Control**
- **Intensity Control**
- **Volume Control**
- **Carries the Battery**

AUDIO SCAN BASIC UNIT

EMDR EQUIPMENT EUROPE

► **Comprises of:**

The Basic Unit and the Headphones plus Tactile Unit



AUDIO SCAN BASIC UNIT WITH TACTILE UNIT

EMDR EQUIPMENT EUROPE

► **Comprises of:**

The Basic Unit and the Headphones plus Tactile Unit plus Protection Bag



SET ONE

AUDIO SCAN BASIC UNIT WITH TACTILE UNIT



Teenager Tactile Units look a bit more interesting

Children's Tactile Units are the Teenager Tactile Units covered with a Rubber duck or similar



TACTILE UNIT TEENAGER/CHILDREN

► **Comprises of:**

Eye Scan Unit with Control box



The Eye Scan;

- Three setting
 - Setting one looping light
 - Setting two looping light with a little stop
 - Setting three end lights are blinking
- Speed adjustable
- Four Colours: Red, Blue, Green, White
- Dimmable

EYE SCAN UNIT

EMDR EQUIPMENT EUROPE

► **Comprises of:**

Eye Scan Unit with Control box, Tripod and Protection bag



SET TWO,

EYE SCAN UNIT WITH TRIPOD AND PROTECTION BAG

EMDR EQUIPMENT EUROPE



Complete Set

- Audio Scan Basic Unit with Headphones
- Tactile Unit
- Special Protection Bag
- Eye Scan Unit
- Tripod
- Special Protection Bag

COMPLETE SET

EMDR EQUIPMENT EUROPE



We recommend to use re-chargeable batteries

This will help the environment and is cheaper to operate

ACCESSORIES, GO GREEN

EMDR EQUIPMENT EUROPE

- ▶ EMDR Equipment is a essential tool for EMDR Therapist
- ▶ All our products are CE Marked
- ▶ All products are covered with a 2 Year Warranty

SUMMERY

EMDR EQUIPMENT EUROPE

▶ **We supply to:**

- ▶ **NHS, MoD, Police Force, Private Hospitals, Charities, Private Practitioners in the UK**
- ▶ **European and world wide to the equivalent organizations**

▶ **All our products are CE Marked**

▶ **All products are covered with a 2 Year Warranty**

SUMMERY

▶ **Other Resources on our website**

- ▶ Books
- ▶ DVD's
- ▶ CD's

▶ **We also support the work of HAP UK & Ireland**

- ▶ Books, Protocols, Scripts, all the income for the HAP Sales is going directly to HAP as a Fundraising for their Projects

▶ **Give Humanity a Chance**

- ▶ Is a CIC to support basic help for Human in need
- ▶ Currently focusing on the help for refugees
- ▶ www.givehumanityachance.org

SUMMERY

EMDR EQUIPMENT EUROPE

more information on our website:
www.emdrequipmenteurope.com



Thank you very much for
listening,

Ulf Jarisch
Managing Director

EMDR Equipment Europe
The Studio
West Wing
Stardens
Tewkesbury Road
Newent
GL18 1LG
01531 820511

EMDR EQUIPMENT EUROPE